



# SuBiTo

# THE SEQUENCING

Typing to the highest resolution

Ready Gene

**SuBiTo**

FluoGene

Software

PCR & Electrophoresis


Ready Plate

DNA-Extraction

Life Science

## SAME SIMPLE AND UNIFORM SET UP FOR ALL HLA-SuBiTo KITS:


- To reduce ambiguities 2 amplifications are carried out in parallel in order to separate many heterozygous allele combinations in the first step
- Gel Electrophoresis
- Purification of the PCR products using ExoSAP-IT™
- 6 sequence reactions per HLA Class I Locus sequence
- 3-5 sequence reactions per HLA Class II Locus sequence
- The combination of forward and reverse reactions takes place according to positive PCR amplification
- Precipitation with IT-Prec-Buffer
- Loading of the capillary sequencer
- Analysis with HLA-HiType software



## NEW KIT IMPROVEMENTS SAVE YOU 3 HOURS ASSAY TIME:

- New chemistry on G1/G2 mixes enable shorter PCR profile for first PCR - takes now only 1 h 45 min
- Shorter protocol for Cycle Sequencing PCR takes now only 1 h 10 min
- Shorter precipitation protocol saves 30 minutes


## CLASS I STRATEGY:



The strategy of 2 group specific amplifications is the ideal solution for high resolution as well as for high throughput. The amplifications are followed by 6 sequencing reactions for exon 2, 3 and 4.

Parallel amplification of exons 5 – 8 as an internal control provides the option of further sequence analysis.

## CLASS II STRATEGY:



The group specific amplifications are followed by 2 sequencing reactions of exon 2 and 1 sequencing reaction of exon 3 .

Parallel amplification of exon 3 or exon 4 respectively as an internal control provides the option of further sequence analysis.

## COMPLETE KIT CONTENTS:

- All components for amplification including IT-Taq Polymerase
- Exo-SAP-IT™ for purification of the PCR products
- Gene locus and exon-specific sequencing mixes incl. BigDye® chemicals
- IT-Prec Buffer for precipitation of the sequencing products
- 2 Second Level Sequencing Mixes for resolution of common ambiguities are directly included in the kit


## SECOND LEVEL SEQUENCING:

Due to the HLA-SuBiTo strategy of separating HLA allele groups by two initial amplifications, common cis-trans ambiguities are strongly reduced. However, in case the groups could not be separated, possible ambiguities can be resolved with the ready-to-use Second Level Sequencing Mixes using the existing PCR product of the standard kit.

There are a total of 31 different Second Level Sequencing Mixes available for HLA class I (10 for HLA-A, 11 for HLA-B and 12 for HLA-C) (some mixes are useable for several loci). For HLA-DRB1 we offer 9 Second Level Sequencing Mixes, 8 mixes for HLA-DQB1 and 4 mixes for HLA-DPB1. The powerful HLA-HiType Software proposes in the case of ambiguities the ideal Second Level Sequencing Mix for resolution.

## ANALYSIS WITH HLA-HiType SOFTWARE:

- Integration of new nomenclature and automated translation in P- and G-groups as well as new NMDP-Codes
- Fully automated data management, no exhausting selection of \*.ab1 Files necessary
- Automated recognition and dynamic display of hemi- and heterozygous sequences
- Several manual and automated editing possibilities
- Diverse alignment options, combined control of electropherogram and alignment
- High safety due to a defined 4 eyes principle (technical and medical validation)
- Complete history with a traceability of each step



## HLA-SuBiTo KITS ARE AVAILABLE FOR 25 TYPINGS EACH:

ARTICLE NO.	HIGH RESOLUTION KITS	ARTICLE NO.	HIGH RESOLUTION KITS
002 090 025*	HLA-SuBiTo A high res	002 093 025*	HLA-SuBiTo DRB1 high res
002 091 025*	HLA-SuBiTo B high res	002 094 025*	HLA-SuBiTo DQB1 high res
002 092 025*	HLA-SuBiTo C high res	002 095 025*	HLA-SuBiTo DPB1 high res
002 096 025*	HLA-SuBiTo ABCDRDQ high res		

\*All articles are also available without BDTs under the same article no. code with /1.

002 090 0SW	HLA HiType Software main licence
002 090 0AL	HLA HiType additional workstation licence



Niederhöchstädter Straße 62  
D-61476 Kronberg/Taunus Germany

Tel. +49 (0)6173- 6079- 30  
Fax: +49 (0)6173- 6079-50

E-mail: info@inno-train.de  
Web: www.inno-train.de