

3DHISTECH


THE DIGITAL PATHOLOGY COMPANY

TRANSITION TO DIGITAL PATHOLOGY

Pannoramic Software 2.0

Increased convenience of use with the new control software for Pannoramic scanners

The new control software for the Pannoramic scanner range from 3DHISTECH, offering live teleconsultation and increased convenience thanks to additional features.



*Convenience
Live teleconsultation*

Key Features

New software architecture:

- Microscope module enabling live image viewing (option to embed third-party image plugins)
- Flexible switching between brightfield and fluorescent mode
- Scanning Planner module - for several scanning profiles in walk-away operation

Remote access for multiple users to enable live group teleconsultation

Pre-assignment of several scanning profiles to each slide (up to the full capacity of the scanner) - including mixed fluorescence and brightfield scan within one scanning batch process


Separate scanning option for tissue microarrays

Customizable user interface themes

Prescan-preview function for all slides

One-click restoration of any previous scanner settings

Log file attached to automated e-mail notification in case of any scan event



Production and
development by



3DHISTECH

WWW.3DHISTECH.COM

3DHISTECH LTD.

H-1141 Budapest, Öv u. 3., Hungary

E-mail for general enquiries: info@3dhistech.com

Phone: +36-1-211-2130

Fax: +36-1-211-2129